

For further inquiries, please contact  
Smile Signature

Smile Signature at Ratchada  
**02 693 6933**

Smile Signature at Phahonyothin  
**02 930 1812**

Smile Signature at Siam Square  
**083 095 0219**

Smile Signature at Ratchayothin  
**080 880 1700**

Smile Signature at Srinagarindra  
**02 721 9102**

Smile Signature at Central West Gate  
**02 194 2733**

Smile Signature at Ramintra  
**02 943 2955**

Smile Signature Onnut  
**02 081 8377**

Smile Signature Rangsit  
**02 958 5600**



Post-Operative  
Instructions for  
Invisalign® aligners



Follow me:

-  SmileSignature
-  SmileSignature
-  SmileSignature

QR CODE  




[www.facebook.com/SmileSignature](http://www.facebook.com/SmileSignature)  
[www.SmileSignature.com](http://www.SmileSignature.com)



## Post-Operative Instructions for Wearing and Using Your invisalign® Aligners



1. Wear your aligners per your doctor's instructions, usually 20–22 hours per day.
2. Wash your hands thoroughly with soap and water before handling your aligners.
3. Handle only ONE aligner at a time.
4. Rinse your aligners when removing them from the packaging.

### Aligner Insertion

1. Make sure you have the proper aligner—the upper for your top teeth and the lower for your bottom teeth.
2. To help avoid confusion, each aligner is engraved with your unique case number, a “U” for upper and an “L” for lower, followed by the stage number.
3. You may insert either the upper or lower aligner first. When inserting each aligner, gently push the aligners over your front teeth. Then apply equal pressure, using your fingertips, to the tops of your left and right molars (back teeth) until the aligner snaps into place.
4. Do NOT bite your aligners into position. This may damage them.

**Note:** If you experience sharp pain or significant discomfort, discontinue use of the aligners and contact your doctor.

### Aligner Removal

1. On one side of your mouth, use your fingertip on the inside of your back molar to slowly pull the aligner from your molars.
2. Repeat this process on the other side of your mouth before trying to completely remove the aligner.
3. Once aligner is disengaged from the molars on both sides of your mouth, you should be able to slowly work your way forward gently prying the aligner away from your teeth with your fingertips.

### Note:

- Immediately rinse aligner with water, shake off excess water, and store your aligners in the protective case provided with your starter kit.
- To help prevent damage to the aligners, avoid unnecessary removal.
- Take care in removing your aligners, especially if multiple attachments are being used.
- Do not use excessive force to bend or twist an aligner to get it off.
- DO NOT use any sharp object to remove your aligners.
- **Consult with your doctor if your aligners are extremely difficult to remove.**

### Daily Care and Maintenance of Your invisalign® aligners

1. Clean your aligners prior to each insertion. Use a soft bristle toothbrush with water and a small amount of toothpaste. You may find it easiest to clean the outside of your aligners by brushing them while they are still on your teeth, then remove your aligners to clean the inside surfaces.

**Note:** Be sure to rinse each aligner thoroughly with water after each cleaning. We also recommend that you use the Invisalign aligner cleaning products once a week, or as needed.

2. Do NOT use denture cleaners to clean aligners. Do NOT soak them in mouthwash. These products can damage the surface of the aligner, causing it to become dull and more visible. (Pic)

### Proper Oral Hygiene

1. Remove your aligners for eating and drinking. (You do not need to remove your aligners to drink cool water.)
2. Brush and floss your teeth after each meal or snack prior to reinserting your aligners. If you don't have access to your cleaning system or a toothbrush, you can simply rinse your mouth, and then clean your aligners by holding them under warm running water. It's not the best way to clean, but it works in a pinch. Be sure to thoroughly clean your aligners at your earliest convenience.
3. If you have any questions regarding hygiene techniques, please consult your doctor.
4. Regular dental checkups and cleaning are recommended for the continued health of your teeth and gums. Storing your Invisalign Aligners.

We recommend that you store your Invisalign aligners in a case when they are not in your mouth. This will help protect them from loss and damage.

**Note :** Most people experience tooth soreness for a few days after starting each new stage. This is normal. It is a sign that the Invisalign aligners are working, moving your teeth to their final destination. This soreness should gradually go away a couple of days after inserting the new aligner in the series. If it doesn't, promptly call your doctor.